

Accessories List



Export

Cryo 6 Skincooling unit



Accessories Cryo 6

95.853.411 Neo light treatment tube Cryo 6, 2 m
95.853.611 Neo light treatment tube Cryo 6, 3 m



95.852.210 Treatment tube light Cryo 6, woven version, 2 m
95.854.410 Treatment tube light Cryo 6, woven version, 3 m




93.852.630 Arm for light treatment tube Cryo 6



65.853.112 Accessory storage dish Cryo 6



65.851.610 Drip tray (condensed water)



CryoMini Skincooling unit4009 **CryoMini 230V/50-60Hz**4011 **CryoMini 120V/50-60Hz**4004 **CryoMini 100V/50-60Hz****Accessories CryoMini**

95.852.210 Treatment tube light, woven version, 2 m

95.854.410 Treatment tube light, woven version, 3 m



93.852.660 Arm for light treatment tube CryoMini



80.026.510 Drip tray (condensed water)



93.493.220 Plate Solta Clear&Brillante for CryoMini



Accessories Cryo 6 / CryoMini / Cryo 5

95.372.620	Adapter for laser handpiece with clip (ø until 20 mm)	
65.373.510	Adapter for laser handpiece without clip	
65.372.631	Adapter flat for laser handpiece without clip	
95.372.710	Clip (ø until 20 mm) for laser handpiece without adapter	
95.373.110	Clip (ø 22-28 mm) for laser handpiece without adapter	
95.373.310	Clip (ø 30-38 mm) for laser handpiece without adapter	
63.373.110	Adapter LightSheer (LightSheer handpiece not included)	
95.855.110	Adapter Mydon	
65.855.010	Adapter Arion	
95.855.610	New hook for Arm Cryo 6 / CryoMini	
87.504.010	Adapter Candela GMAX	
87.504.020	Adapter Candela MGLase/GYAG	

Accessories Cryo 6 / CryoMini / Cryo 5	
95.443.010	Hybrid Air Clip for Lutronic "Clarity" <i>(Clarity handpiece not included)</i>
95.853.620	Lumenis Tube for UltraPulse, 3 meters <i>(Ultrapulse handpiece not included)</i>



Accessories Cryo 5	
95.372.210	Treatment tube light
95.372.212	Treatment tube light special length 3m
65.370.620	Holder for light treatment tube



ZWave Shockwave for Aesthetics treatments including SysCartPro	
5530-X01	ZWavePro Device with accessories:
5530-x02	ZwavePro Device with 1 Handpiece <i>Warranty for the shockwave generator of ZWave: 2 Million shocks</i> <i>Applicator heads of ZWave are consumables and are not included in the warranty</i>



ZWave Shockwave table unit for Aesthetics treatments without SysCartPro	
5530-X03	ZWavePro Device table unit with accessories: 2x Handpiece 2x Holder for Handpiece 3x Applicator head (39 mm) 1x Lotion for treatment 1x Footswitch 10x Silicone caps 1x Template
5530-X04	ZWavePro Device table unit with 1 Handpiece <i>Warranty for the shockwave generator of ZWave: 2 Million shocks</i> <i>Applicator heads of ZWave are consumables and are not included in the warranty</i>



ZWave Shockwave for Aesthetics treatments

5518-X01

ZWave Device with accessories:

- 1x Handpiece
- 1x Holder for Handpiece
- 3x Applicator head (39 mm)
- 1x Lotion for treatment
- 1x Footswitch
- 1x Transport case
- 1x Template
- 20x Silicone caps

*Warranty for the shockwave generator of ZWave: 2 Million shocks
Applicator heads of ZWave are consumables and are not included in the warranty*



ZWave with Move Cart

5518-X02

ZWave Device with accessories:

- 1x Handpiece without applicator head
- 1x Holder for Handpiece
- 3x Applicator head (39 mm)
- 1x Lotion for treatment
- 1x Footswitch
- 1x Transport case
- 1x Template
- 20x Silicone caps

*Warranty for the shockwave generator of ZWave: 2 Million shocks
Applicator heads of ZWave are consumables and are not included in the warranty*



914

Move Cart

Total Package

SysCartPro

9160

for ZWavePro



Cart Move

914

Move Cart
Including drawer and glass shelf



96.116.010 Glas plate for move cart (optional)

Cart Move SL

915

Move SL
Including glass shelf



Accessories ZWavePro		
5511	Handpiece ZWave (Silicone cap not included)	
5413	Handpiece enPuls Version 2.0	
94.130.411	Foot pedal switch with cable ZWave	
93.133.530	Applicator head 39 mm	
50.500.038	Lotion for application	
65.137.110	Silicone cap (1pc) for applicator head 39 mm	
93.130.411	Holder for handpiece	
95.800.610	Swivelling base for setting on table	

Accessories ZWave		
5511	Handpiece ZWave Version 2.0 (Silicone cap not included)	
94.130.410	Foot pedal switch with cable ZWave	
93.133.530	Applicator head 39 mm	
50.500.038	Lotion for application	
87.053.010	Transport case	
65.137.110	Silicone cap (1pc)	
93.130.411	Holder for handpiece	
65.800.410	Touch pen	

Hyaluronic Filler - ZFill

50.900.130 ZFill refresh²

Product	Hyaluronic acid gel with glycerin
Treatment Option	Mesotherapy, lasting hydration and revitalisation
Areas to be treated	Face, neck, neckline, backs of the hands, upper arms
Hyaluronic acid concentration	18 mg/ml
Stabilizer / crosslinker	Glycerin
Degree of cross-linking	-
Molecular weight	2,5 mio. daltons
Injection depth	Upper dermis
Duration of effect	3 – 4 months after
	3 treatments in the first two months
Packing	1 x 1 ml
Needle	2 x 30G ½

50.900.110 ZFill contour²

Product	Hyaluronic acid-gel-implant
Treatment Option	Correction of superficial and deeper facial lines, lip contours and lip volume
Areas to be treated	Brow lines, glabella lines, periorbital wrinkles, nasolabial lines, periorbital rhytides, marionette lines, lip volume and lip contours
Hyaluronic acid concentration	23 mg/ml
Stabilizer / crosslinker	BDDE
Degree of cross-linking	+
Molecular weight	2,5 mio. daltons
Injection depth	Middle dermis
Duration of effect	6 – 9 months
Packing	1 x 1 ml
Needle	2 x 27G ½

50.900.120 ZFill deep²

Product	Hyaluronic acid-gel-implant
Treatment Option	Correction of deep lines and wrinkles, sunken areas of skin, facial contours and increased volume
Areas to be treated	Cheeks, nasolabial lines, marionette lines, chin, sunken facial areas
Hyaluronic acid concentration	23 mg/ml
Stabilizer / crosslinker	BDDE
Degree of cross-linking	++
Molecular weight	3 mio. daltons
Injection depth	Lower dermis
Duration of effect	8 – 12 months
Packing	1 x 1 ml
Needle	2 x 27G ½



ZCryo Marketing documents

10.101.854 Overview brochure Aesthetic division english



10.101.863 Brochure ZCryo Cryo 6 Laser english



10.101.862 Brochure ZCryo CryoMini english
(max. 50 pcs)



10.101.981 Cryo CD Version 2 (studies, pictures...)
(max. 10pcs)



ZWavePro + ZWave Marketing documents

10.101.815	Patient brochure english (max. 50 pcs)		
10.101.796	Studies brochure english + german (max. 30 pcs)		
10.101.810	DVD-Video german and english (max. 10 pcs)		

ZWavePro Marketing documents

10.102.178	Brochure ZWave ^{Pro} english (max. 50 pcs)		
------------	--	--	--

ZWave Marketing documents

10.101.378	Brochure ZWave english (max. 50 pcs)		
------------	---	--	---

Aesthetic Division Marketing Materials

10.103.075	USB Key 8 GB		
10.103.057	Zimmer Aesthetic Division Notice bloc with pen		
10.103.085	Zimmer Aesthetic Division Cup		
10.103.049	Zimmer Aesthetic Division Bag		
10.103.084	Zimmer Aesthetic Division Umbrella		
10.103.067	Infrared Thermometer with Aesthetic logo		

ZFill Marketing documents

10.101.990 Patient brochure



10.101.970 Doctor brochure



10.101.385 DVD-Video „the new virtually painless concept“



10.103.102 Roll-up for waiting-room english



10.101.384 Poster DIN A2



10.102.136 Patient questionnaire



10.102.157 Product Cards ZFill refresh²



10.102.159 Product Cards ZFill contour²



10.102.158 Product Cards ZFill deep²



Delivery Terms

- Payment in advance
- Delivery ex works
- Units delivered with a 2 years warranty on parts
- During the warranty time we deliver the necessary spare parts at no costs for repair to the distributor. The dealer is responsible for repair and servicing
- Please note the ZFill handling and storage instructions!

Zimmer MedizinSysteme GmbH
Junkersstraße 9
D-89231 Neu-Ulm, Germany
Tel. +49 7 31. 97 61-291
Fax +49 7 31. 97 61-299
export@zimmer.de
www.zimmer-aesthetics.com